



12555 Loma Rica Drive, Grass Valley, CA 95945
 530-273-8888 • 800-WE PIPET • FAX 530-273-2586

Bottle-Top Dispensers

RETAIL PRICE LIST
 Effective: January 1, 2008

CAT. NO.	DESCRIPTION	Price
<p>StatMatic Dispenser StatMatic push-button dispensers (bottle-top). Volume range: 50µl to 2.5ml. Polypropylene and glass construction, with Teflon® plungers. Each instrument includes: 1 each 20, 24 and 28mm threaded adapters; 7" inlet tube; 4 ½" outlet delivery tube; 25" remote delivery tube; volume calibration rings and instructions.</p>		
401	StatMatic – volumes: 50, 100, 200 and 250µl	\$175.50
402	StatMatic – volumes: 250, 300, 400 and 500µl	\$175.50
981	StatMatic – volumes: 0.5, 0.8 1.0, 2.0 and 2.5ml	\$175.50


StatMatic I 8 Channel Microplate Dispenser (bottle top): Polypropylene and glass construction with Teflon® plunger. Each instrument includes: 3 each dispense volume spacers, 1 each 20mm, 24mm and 28mm threaded adapters, 1 each flexible remote outlet tubing, 1 each long inlet tubing, 1 each 8-port manifold, manifold vertical holder and instructions. Volumes dispensed: .8, 1.6, 2.0 or 2.4ml. With manifold: 100, 200, 250 or 300µl per well of a 96 well plate. Manifold is 8-port spacing to fit a 96 well microplate.


982	StatMatic I 8 Channel Microplate Dispenser	\$210.40
-----	--	----------

StatMatic II 8 Channel Microplate Dispenser: Polypropylene and glass construction with Teflon plunger. Each instrument includes: 500mL plastic bottle with 28mm cap assembly which includes bottle inlet tubing and luer fitting, 3 dispense volume spacers, one 30" long Tygon remote inlet tubing and instruction manual. Volumes dispensed: .8, 1.6, 2.0 or 2.4mL. With manifold: 100, 200, 250 or 300µL per well of a 96 well plate. Manifold is 8-port spacing to fit a 96 well microplate.


980	StatMatic II 8 Channel Microplate Dispenser	\$251.50
-----	---	----------

See reverse side for Millimatic Dispensers and
 Accessories for all bottle-top dispensers





CAT. NO.	DESCRIPTION	Price
Millimatic Dispenser		
158	Millimatic Dispenser, inert fluorocarbon and glass construction. Range: 1 to 14.9ml Glass Plunger – 53mm cap and 28mm adapter	\$222.25
149	Millimatic Dispenser, inert fluorocarbon and glass construction. Range: 0.25 to 3.0ml Glass Plunger – 53mm cap and 28mm adapter	\$222.25
803	Millimatic 8-Channel Microplate Dispenser (bottle top) Range: 0.25 to 3mL Inert fluorocarbon, Polypropylene and glass construction. Each instrument includes: glass plunger, 53mm jar cap and 28mm bottle adapter; bottle/jar inlet tubing; one 16" long flexible remote outlet tubing; one 8-port manifold with fittings and adapters; manifold holder and instruction insert.	\$239.50
Accessories		
700	StatMatic Repair Kit: Inlet and outlet tubing, remote tubing, and ceramic ball.	\$69.70
699	Millimatic Repair Kit: Inlet and outlet tubing, remote tubing, ceramic ball, outlet seal cap, calibration scales with knob screw, 28, 32 and 38mm adapters	\$123.60
109	Millimatic 950cc amber glass bottle with 53mm threaded neck	\$37.30
156	Millimatic 500cc amber glass bottle with 53mm threaded neck	\$63.30
790	125cc HDPE bottle with a 3.65" tubing inlet and an adapter to fit StatMatic and StatMatic I Plate Washer/Filler	\$63.30
822	Manifolds, 8-port, non-sterile/disposable, P.P., 25/box	\$119.25
925	Fixed outlet tubing for StatMatic, 5/pkg	\$19.00



Note: Bottles are not included with dispenser.
To select appropriate adapter, measure the outside diameter (O.D.) of the bottle threads in millimeters.
All products are made in the U.S. A.
Warranty: One year, parts and labor on dispensers.