StatPet

FEATURES

- > Convenient, accurate and cost effective
- > Single volume air displacement pipettes
- > Custom colors with your name and logo

CUSTOM OPTIONS

- > Available for OEM customers only
 - Dispense Volume: 20 μ L 1 mL

Unit Color

Handle and Knob Printing

Packaging

OEM REQUIREMENTS

> A volume of 100 units allows custom options




BASE MODEL

This product is available as OEM only



(Shown with optional tip)



SPECIFIC/	ATIONS													
Model	20 µL	50 μL	100 μL	150 μL	200 μL	250 μ	L	300 μL	500 μL	700	μL	800 μL	1 mL	
Tips	200 μL									1mL				
Liquid Co	ntact Mate	rials :	Polypropylene	(disposable	tip)									
FUNCTIO	NAL PARAN	/IETERS												
Precision (full stroke) 2% CV or better			er											
Accuracy (full stroke) ±3%														
200			Volume	Volume		± % Accuracy		% Precision (CV)						
			200 μL - 1 mL	•		3		2						
			50 μL - 199 μ 20 μL - 49 μL		5 7		2 3							
Weight .72 3.1			.72 oz [22.24	95 in. [2.02 cm] 2 oz [22.24 gm] (instrument only) 1 oz [88 gm] (shipping weight)										
	OPTIONS													
Unit Colo	or (s) - Body,	Knob and C	Сар											
Tan (<i>PMS# 4535C</i>) Gray (<i>PMS# 4</i>		MS# 413U)	Teal (PMS# 3145C)		Rose (<i>PMS# 222</i>		22U)	Blue (PMS#	285C)	Ye	llow (PMS# 2	129C)		
Barrel Co	lor(s)													
White Bluish—White														
PAD PRINTING OPTIONS - HANDLE				HOT ST	HOT STAMPING OPTIONS - KNOB									
White or Black Ink Color				White c	White or Black Ink Color									
45 in. [11.4 mm]														
	.40	Print Area in. x 1.89 in mm x 48 mn		Print Area ∅.30 in. [Ø 7.62 mm]										

One time set up charge may apply.

MINIMUM OEM ORDER / SHIPMENT REQUIREMENTS

100 unit minimum order

25 unit minimum per shipment

PRECISION DISPENSING MICROPIPETTE FOR OEMS WORLDWIDE

TriContinent by Gardner Denver

tricontinent.com



12555 Loma Rica Drive Grass Valley, CA 95945 USA Tel: 800-937-4738 Fax: 530-273-2586 liquidhandling.tcs@gardnerdenver.com Gardner Denver Thomas GmbH Livry-Gargan-Str. 10 82256 Fürstenfeldbruck Germany Tel: +49 8141 2280 0 Fax: +49 8141 8892136 thomas.de@gardnerdenver.com



by Gardner Denver

Gardner Denver Hong Kong, Limited Unit 1317-1318 Delta House, 3 On Yiu Street Siu Lek Yuen, Shatin, New Territories Hong Kong Tel: +852 26903502 Fax: +852 27924598 thomas.hk@gardnerdenver.com

The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as an aid to help in the selection of TriContinent products. It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability in connection there with. TriContinent does not warrant, guarantee or assume any obligation or liability in connection with this information.

Photos of products pictured in this catalog do not necessarily represent a specific model number. To obtain further information for custom options, contact your local TriContinent office.

Printed in USA Form No. MKT90012 H © Gardner Denver TriContinent. All rights reserved.