

C-Series Rotary Valve

FEATURES

- > 2-year warranty against defects in materials and workmanship
- > Long-life
- > Exceptional Value
- > Available in a variety of inert materials
- > No solenoid valve limitations (leaking, clogging, low CVs)
- > No ceramic valve limitations (seizing)

CONFIGURED OPTIONS

- > Materials
 - Body PEEK
PCTFE
 - Rotary Plug UHMWPE
PTFE
- > Valve Ports
 - "Y" 120° 3-port
 - 90° 4-port
 - Distribution 4-port
 - T 3-port
 - T 4-port
 - Distribution 6-way

OEM OPTIONS




- > Custom configurations possible

BASE MODEL

Compatible with C-Series Pumps and Rotary Valve Controllers



SPECIFICATIONS

Description	Materials (Body/Plug)	Catalog Number	Flow Path (Blue)
"Y" 120° 3-port	PEEK/UHMWPE	10020	
	PEEK/PTFE	10021	
	PCTFE/UHMWPE	10022	
	PCTFE/PTFE	10023	
90° 4-port	PEEK/UHMWPE	10040	
	PEEK/PTFE	10041	
	PCTFE/UHMWPE	10042	
	PCTFE/PTFE	10043	
Distribution 4-Port	PCTFE/UHMWPE	10062	
	PCTFE/PTFE	10063	
T 3-port	PEEK/UHMWPE	10115	
T 4-port	PEEK/UHMWPE	10125	
Distribution 6-way	PCTFE/UHMWPE	10106	

ROTARY VALVES FOR SYRINGE PUMPS

tricontinent.com



TriContinent Scientific, Inc.
 12555 Loma Rica Drive
 Grass Valley, CA 95945
 USA
 Tel: 800-937-4738
 Fax: 530-273-2586
liquidhandling.tcs@gardnerdenver.com

Gardner Denver Thomas GmbH
 Livry-Gargan-Str. 10
 82256 Fürstentfeldbruck
 Germany
 Tel: +49 8141 2280 0
 Fax: +49 8141 8892136
thomas.de@gardnerdenver.com

Gardner Denver Hong Kong, Limited
 Unit 1317-1318 Delta House, 3 On Yiu Street
 Siu Lek Yuen, Shatin, New Territories
 Hong Kong
 Tel: +852 26903502
 Fax: +852 27924598
thomas.hk@gardnerdenver.com



The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as an aid to help in the selection of TriContinent products. It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability in connection there with. TriContinent does not warrant, guarantee or assume any obligation or liability in connection with this information.

Photos of products pictured in this catalog do not necessarily represent a specific model number. To obtain further information for custom options, contact your local TriContinent office.

Printed in USA Form No. MKT90034 C © Gardner Denver TriContinent. All rights reserved.