Prime Series Syringes



FEATURES

- > Seal materials include the highest quality PTFE or UHMWPE
- > Compatible with industry standard 3 cm and 6 cm stroke syringe pumps
- > Highest quality borosilicate 3.3 glass barrels
- > Bores held to extremely tight tolerances

CONFIGURED OPTIONS

Seal Materials
UHMWPE - for best abrasion resistance and longest life
PTFE - for best chemical resistance

OEM OPTIONS

> Custom volumes and configurations available



SPECIFICATIONS

Volume	50 μL	100 μL	250 μL	500 μL	1.0 mL	2.5 mL	5.0 mL	12.5 mL*
Seal Materials	UHMWPE	(black) or P	ΓFE (white)		_			
Barrel Material	Borosilicat	e 3.3 Glass			•			
Fitting Interface	1/4-28				•			
Stroke	3 cm							

^{*}PTFE only



Volume	50 μL	100 μL	250 μL	500 μL	1.0 mL	2.5 mL	5.0 mL	
Seal Materials	UHMWPE (black) or PTFE (white)							
Barrel Material	Borosilicat	e 3.3 Glass						
Fitting Interface	1/4-28							
Stroke	3 cm							



Volume	50 μL	100 μL	250 μL	500 μL	1.0 mL	2.5 mL	5.0 mL	10.0 mL	25.0 mL
Seal Materials	UHMWPE	(black) or P	TFE (white)						
Barrel Material	Borosilicate 3.3 Glass								
Fitting Interface	1/4—28				•				
Stroke	6 cm				•				





by Gardner Denver

tricontinent.com



TriContinent Scientific, Inc. 12740 Earhart Ave

Auburn, CA 95602 USA Tel: 800-937-4738 Fax: 530-273-2586

liquidhandling.tcs@gardnerdenver.com

Gardner Denver Thomas GmbH

Livry-Gargan-Str. 10 82256 Fürstenfeldbruck Germany Tel: +49 8141 2280 0 Fax: +49 8141 8892136

thomas.de@gardnerdenver.com

Gardner Denver Hong Kong, Limited

Unit 1317-1318 Delta House, 3 On Yiu Street Siu Lek Yuen, Shatin, New Territories Hong Kong Tel: +852 26903502

Fax: +852 27924598

thomas.hk@gardnerdenver.com

The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as an aid to help in the selection of TriContinent products. It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability in connection there with. TriContinent $does \ not \ warrant, \ guarantee \ or \ assume \ any \ obligation \ or \ liability \ in \ connection \ with \ this \ information.$

Photos of products pictured in this catalog do not necessarily represent a specific model number. To obtain further information for custom options, contact your local TriContinent office.